

Southern Nevada Public Health Laboratory
TEST MENU

Complete Blood Count (CBC)

Method	Sysmex® XN-2000 multi-parameter automated hematology analyzer - the hydrodynamic focusing (DC Detection), flow cytometry method (semiconductor laser), and SLS-hemoglobin method.																																																																										
CPT Code	85027																																																																										
Specimen	Whole blood – Required Volume: 1mL RBT Micro tube – 250uL Minimum specimen volume: 1mL																																																																										
Collect in	EDTA-2K or EDTA-3K anticoagulant tube																																																																										
Labeling	Label with patient's full name, date of birth, collection date, and time.																																																																										
Processing	Processing Serum: Allow refrigerated samples to come to room temperature and mix well before analysis. Do not place CBC and Diff samples on a mechanical rocker. Constant rocking may alter white cell membranes, resulting in false interpretive messages.																																																																										
Transport	Store and transport refrigerated (2-8°C) within 48 hours Room temperature stability is 24 hours. Keep tube tightly capped for storage.																																																																										
Sample Rejection	<ul style="list-style-type: none"> - Mislabelled/unlabeled specimen - Specimen quantity not sufficient for testing (QNS) - Specimen contains fibrin or clots. - Excessive platelet clumping - Substandard mixing or collection 																																																																										
Requisition	SNPHL Request Forms																																																																										
Turn Around Time	24-48 Hours																																																																										
Results	<p>Reference Range:</p> <p>Adult:</p> <table border="1"> <thead> <tr> <th>Parameter</th><th>Males (21-99)</th><th>Females (21-99)</th><th>Units</th></tr> </thead> <tbody> <tr> <td>WBC</td><td>3.18-12.74</td><td>2.74-10.9</td><td>x103/mcL</td></tr> <tr> <td>RBC</td><td>4.08-5.80</td><td>3.62-5.76</td><td>X106/mcL</td></tr> <tr> <td>HGB</td><td>12.3-17.2</td><td>10.7-16.7</td><td>gm/dL</td></tr> <tr> <td>HCT</td><td>37-52</td><td>33-50</td><td>%</td></tr> <tr> <td>MCV</td><td>83.0-97.0</td><td>79.4-99.5</td><td>fL</td></tr> <tr> <td>MCH</td><td>27.3-32.4</td><td>24.8-33.8</td><td>Pg</td></tr> <tr> <td>MCHC</td><td>31.3-35.1</td><td>30.6-34.9</td><td>gm/dL</td></tr> <tr> <td>RDW-CV</td><td>11.4-14.7</td><td>11.0-15.1</td><td>%</td></tr> <tr> <td>RDW-SD</td><td>35.3-50.7</td><td>37.6-47.5</td><td>fL</td></tr> <tr> <td>IPF</td><td>1.23-8.70</td><td>Same as Males</td><td>%</td></tr> <tr> <td>PLT</td><td>140-400</td><td>Same as Males</td><td>x103/mcL</td></tr> <tr> <td>MPV</td><td>7.2-12.6</td><td>Same as Males</td><td>fL</td></tr> <tr> <td>NRBC%</td><td>0.0-0.2</td><td>Same as Males</td><td>/100WBC</td></tr> <tr> <td>NRBC#</td><td>0.0-0.01</td><td>Same as Males</td><td>x103/mcL</td></tr> <tr> <td>RETIC %</td><td>0.51-1.81</td><td>0.50-1.70</td><td>%</td></tr> <tr> <td>RETIC ABS</td><td>2.60-9.50</td><td>1.64-7.76</td><td>104/mcL</td></tr> <tr> <td>IRF</td><td>0.0-3.4</td><td>0.0-3.5</td><td>%</td></tr> </tbody> </table>			Parameter	Males (21-99)	Females (21-99)	Units	WBC	3.18-12.74	2.74-10.9	x103/mcL	RBC	4.08-5.80	3.62-5.76	X106/mcL	HGB	12.3-17.2	10.7-16.7	gm/dL	HCT	37-52	33-50	%	MCV	83.0-97.0	79.4-99.5	fL	MCH	27.3-32.4	24.8-33.8	Pg	MCHC	31.3-35.1	30.6-34.9	gm/dL	RDW-CV	11.4-14.7	11.0-15.1	%	RDW-SD	35.3-50.7	37.6-47.5	fL	IPF	1.23-8.70	Same as Males	%	PLT	140-400	Same as Males	x103/mcL	MPV	7.2-12.6	Same as Males	fL	NRBC%	0.0-0.2	Same as Males	/100WBC	NRBC#	0.0-0.01	Same as Males	x103/mcL	RETIC %	0.51-1.81	0.50-1.70	%	RETIC ABS	2.60-9.50	1.64-7.76	104/mcL	IRF	0.0-3.4	0.0-3.5	%
Parameter	Males (21-99)	Females (21-99)	Units																																																																								
WBC	3.18-12.74	2.74-10.9	x103/mcL																																																																								
RBC	4.08-5.80	3.62-5.76	X106/mcL																																																																								
HGB	12.3-17.2	10.7-16.7	gm/dL																																																																								
HCT	37-52	33-50	%																																																																								
MCV	83.0-97.0	79.4-99.5	fL																																																																								
MCH	27.3-32.4	24.8-33.8	Pg																																																																								
MCHC	31.3-35.1	30.6-34.9	gm/dL																																																																								
RDW-CV	11.4-14.7	11.0-15.1	%																																																																								
RDW-SD	35.3-50.7	37.6-47.5	fL																																																																								
IPF	1.23-8.70	Same as Males	%																																																																								
PLT	140-400	Same as Males	x103/mcL																																																																								
MPV	7.2-12.6	Same as Males	fL																																																																								
NRBC%	0.0-0.2	Same as Males	/100WBC																																																																								
NRBC#	0.0-0.01	Same as Males	x103/mcL																																																																								
RETIC %	0.51-1.81	0.50-1.70	%																																																																								
RETIC ABS	2.60-9.50	1.64-7.76	104/mcL																																																																								
IRF	0.0-3.4	0.0-3.5	%																																																																								
Reported	System generated fax; Electronic transmission																																																																										
Limitations	Specimens must be free of clots and fibrin strands. Severely hemolyzed samples (<i>in vitro</i>) falsely decrease RBC and hematocrit.																																																																										
Contact	Southern Nevada Public Health Laboratory – Clinical Laboratory Department																																																																										