Southern Nevada Public Health Laboratory TEST MENU

Alkaline Phosphatase (ALP)

Method	Abbott Alinity c-series Para-nitrophenyl Phosphate
Wicthou	Abbote Amin'ty & Series Fura merophenyi Friosphate
CPT Code	84075
Specimen	Serum – Required Volume: 2-4 mL
•	Minimum specimen volume: 1mL
Collect in	Serum Separator Tube (SST)
Labeling	Label with patient's full name, date of birth, collection date, and time.
Processing	Processing Serum: Allow blood to completely clot (usually within 30 minutes). Serum may be
	centrifuged at the speed and time recommended by the manufacturer.
Transport	Store and transport refrigerated (2-8°C) within 7 days.
Sample Rejection	- Mislabeled/unlabeled specimen
Sumple Rejection	- Specimen quantity not sufficient for testing (QNS)
	opening quantity necessition costing (Q.16)
Requisition	SNPHL Request Forms
Turn Around Time	24-48 Hours
Results	Reference Range:
	Adult Male: (> 20 years) 40 to 150 U/L
	Adult Female: (>15 years) 40 to 150 U/L
Reported	System generated fax; Electronic transmission
1* 1*	
Limitations	
Contact	Southern Nevada Public Health Laboratory – Clinical Laboratory Department
Contact	300 the first vevada i abilic fical tri Laboratory — Cillical Laboratory Department