## Southern Nevada Public Health Laboratory TEST MENU

## **Candida auris Susceptibility**

Method	Broth Microdilution and Etest
CPT Code	
Specimen	Candida auris isolates submitted by local providers. C. auris suspected of causing clinical illness.
Collect in	
Labeling	Label transport medium (slant or plate) with patient's first and last name, date and time of collection
Processing	Indicate specimen source on requisition form. Record patient symptomatology on request form.
Transport	Sabouraud Dextrose Agar Slant/Plate
Sample Rejection	<ul> <li>Mislabeled/unlabeled specimen or Requisition form/isolate discrepancies</li> <li>Contaminated culture medium</li> <li>Broken transport medium</li> <li>Leaking specimen</li> <li>Improperly sealed slant or plate</li> </ul>
Requisition	SNPHL Request Forms Testing is performed only at the request of SNHD
Turn Around Time	1-2 weeks
Results	Minimum Inhibitory Concentration, no interpretation
Reported	System generated fax; Electronic transmission
Note	
Contact	Southern Nevada Public Health Laboratory – Microbiology Department